
KIT-UR-G

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
Dimensioni assieme (mm) Assembly dimensions (mm)

KIT-UR-G

Electric 2-finger parallel gripper with plastic cover and capacitor box. The system can be assembled on the UR3, UR5 or UR10 collaborative robots and the included capacitor box allows for a direct connection to the robot wrist. Standard fingers are also included.

Main features


- A unique design for UR 3, UR 5 and UR 10.
- No cables along the robot arm: direct connection to M8 plug connector of the robot.
- Easy to install and without any configuration (plug & play solution).
- Simulated proximity switch functionality embedded for gripping or end of the stroke detection.

How does it work?

Mount the gripper with the included mechanical flange to the wrist of the robot and install the protective cover and the capacitor box. The system can be electrically connected directly to the M8 plug connector at the wrist of the robot (not cables along the robot arm). The system can be operated via digital input and output to/from the tool interface of the UR robot.

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
**Dimensioni (mm)
Dimensions (mm)**


KIT-UR-V

End Of The Arm Tool (EOAT) for vacuum based pick & place operations. The system can be assembled on the UR 3, UR 5 or UR 10 collaborative robots and the included pneumatic actuator and components allow for a direct connection to the robot wrist and to the pneumatic power source.

Main features


- A unique design for UR 3, UR 5 and UR 10.
- No electric cables along the robot arm: direct connection to M8 plug connector of the robot.
- Easy to install and without any configuration (plug & play solution).
- Embedded vacuum generator, valve and solenoid.

How does it work?

Mount the system with the included mechanical flange to the wrist of the robot and connect the main valve inlet to the pneumatic power source. The system can be electrically connected directly to the M8 plug connector at the wrist of the robot (not cables along the robot arm). The system can be operated via digital output from the tool interface of the UR robot.

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

KIT-UR-J

Parallelogram electric gripper with integrated capacitor box. The gripper can be assembled on all the UR collaborative robots (UR3,UR5,UR10) and the included capacitor box allows for a direct connection to the robot wrist. The user can install his/her own fingers connecting to the existing ones. This gripper is suitable for long stroke and fast opening/closing motion and the gripping force is unaffected by the length of the custom fingers. By installing the associated Gimatic URCap the user can configure the device and a list of work pieces (WP) to be gripped just once and simply use a generic grip/release command independently on the internal/external gripping configuration. The plugin also automatically updates the mass distribution of the payload and the final tool center point (TCP) of the gripper. The URCap provided by Gimatic can handle multiple Gimatic's UR KITs and WP definitions simultaneously in a unique application thus simplifying the development process.

Main features


- A unique design for UR3, UR5 and UR10 and e-Series.
- No cables along the robot arm: direct connection to M8 plug connector of the robot.
- Easy to install and without any configuration (plug & play solution).

Simulated proximity switch functionality embedded for gripping or end of the stroke detection.

How does it work?

Mount the gripper with the included mechanical flange to the wrist of the robot. The system can be electrically connected directly to the M8 plug connector at the wrist of the robot (not cables along the robot arm). The system can be operated via digital input and output to/from the tool interface of the UR robot of using the Gimatic's URCap.

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
KIT-UR-QC

The KIT-UR-QC is a manual tool changer kit for UR collaborative robots. It's mainly composed of two parts: one permanently connected to the robot wrist (QCY90-A) and the other permanently connected to a tool (QC90-B). By operating a mechanical lever, the user can dis/connect the two parts allowing for an easy tool replacement. The system can be assembled on the UR 3, UR 5 or UR 10 collaborative robots and it includes both pneumatic and electric connection components useful to connect power sources between robot and tool. The whole system is a plug & play device with direct connection to the robot wrist and it can be used to transfer either pneumatic pressure and vacuum to the tool. Additional QC90-B spare parts are available upon request to create a more complex system where a robot interacts with several tools.

Main features

- A unique design for UR 3, UR 5 and UR 10.
- No electric cables along the robot arm: direct connection to M8 plug connector of the robot.
- Easy to install and without any configuration (plug & play solution).
- Directly compatible with others Gimatic's UR kits like KIT-UR-G and KIT-UR-V.
- Embedded pneumatic fittings and electric connectors.

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How does it work?

Install the system on the robot's wrist using the included mechanical flange and connect the pneumatic fittings and the electrical connectors as required. The system can be connected directly to the M8 connector on the robot's wrist and can therefore be controlled using the standard digital outputs available.